

**Stabilizor®**

**Peptide extraction kit, Tissue**

**User Manual: LC-MS**

• **Introduction**

Analysis of endogenous peptides in biological tissues is easily compromised by the rapid emergence of *ex vivo* protein degradation fragments. Stabilization of tissue samples in the Stabilizor™ system inactivates the enzymes that cause sample degradation and thereby prevent Protein/peptide changes.

This protocol focuses primarily on extraction of water soluble neuropeptides from brain tissue – molecules particularly prone to rapid degradation. However, the protocol is also valid for other sample types.

• **Materials supplied**

- 12 x Maintainer™ Tissue cards
- 12 x 10kDa Cut-off spin filters
- 24 x Tubes for homogenization/centrifugation
- 3 x 45 ml Extraction fluid
- User Manual: LC-MS

• **Storage**

- Store kit at room temperature
- Use Extraction Fluid within 24 months from date of manufacture (see label)
- After opening, store Extraction Fluid at +4°C-+8°C and use within one month

• **Required equipment (not provided)**

- Stabilizor T1 instrument
- Adjustable pipette (1000µl)
- Analytical balance scale
- Homogenization device
- Centrifuge (14,000 x g)
- 1.5 ml centrifuge tubes

• **Safety information**

Concentrations of hazardous compounds in the Extraction Fluid are below those set by REACH Regulation (1907/2006) as requiring Safety Data Sheets. A Safety Data Sheet can be supplied upon request.

• **Filter limitations**

Filters are non-autoclavable.

Do not centrifuge above 14,000 x g.

Filter membranes contain trace amounts of glycerin. If glycerin will interfere with downstream analysis, pre-wash filters by centrifugation of the assembled filter device using dH<sub>2</sub>O or extraction buffer. After washing, do not let filter membrane dry out completely before use.

• **Sample preparation**

Avoid unnecessary freeze-thaw cycles and variations in sample handling. Stabilized samples are not sensitive to enzymatic activity; however non-enzymatic processes such as oxidation can still affect the samples.

• **Technical assistance**

For questions about this kit as well as the use of the Stabilizor instrument and working with heat-stabilized tissue, visit [www.denator.com](http://www.denator.com) for up to date details of application and technical support in your region.

Products are for research use only

## Protocol

### Sample collection

1. Prepare all equipment and solutions prior to sampling.
2. Extract sample and place centrally in a Maintainor Tissue card.
3. Stabilize immediately in the Stabilizor system using the selected program mode.  
*! Samples can be frozen and stabilized at a later stage if preferred*
4. Store stabilized samples in a freezer if they are not to be analyzed directly.  
*! Stabilized samples can be dissected if required. Avoid pieces smaller than 20 mg as some loss may occur in the spin filter.*
5. Transfer samples to pre-weighed spin filter collection tubes.  
*! Keep frozen samples on dry ice while weighing collected tubes.*
6. Weigh samples and calculate sample weight.

### Homogenization

Samples must be thoroughly homogenized to facilitate peptide solubilisation.

Microtip sonication works well when homogenizing stabilized brain tissue.

1. Add extraction buffer, 5µl/mg of sample.
2. Homogenize e.g microtip sonication 3x4 seconds in an ice bath.
3. Centrifuge extract at top speed for approximately 30 minutes to clarify suspension and pellet cell debris.
4. Collect supernatant.

### Protein concentration measurement

*For semi-quantitative measurements*

Use a standard kit such as the BCA™ Protein Assay Kit (Pierce) to part of each sample. Keep other part frozen during analysis. Add extraction buffer to level out any concentration differences.

### Peptide separation

1. Assemble filter and centrifuge tube – insert white end of filter into the centrifuge tube.
2. If necessary, pre-wash filter by adding 50µl extraction buffer and centrifuge at 14,000 x g, for 15 minutes. Remove filter and put it into a new collection tube.  
*! When using pre-washed tubes, use equal volumes of extract from all samples in step 3.*
3. Carefully decant 100 - 500µl, of extract and pipette into the filter reservoir.  
*! Do not touch filter with the pipette tip.  
! Ensure that particulate matter does not enter the reservoir as this may clog the membrane.*
4. Centrifuge at 14,000 x g for 60 minutes.
5. If extracting from small samples (<20 mg), wash filter once by adding 25µl of extraction buffer and centrifuge again at 14,000 x g for 15 minutes.
6. Remove filter from the centrifuge tube.
7. Store filtrate containing peptide extract for subsequent use.

### LC-MS and analysis

Subsequent downstream analysis is not affected by using the Stabilizor system. Conventional LC-MS protocols can be followed without alteration.

Reference:

*Heat Stabilization of the Tissue Proteome: A New Technology for Improved Proteomics*

*M. Svensson et al. J. Proteome Res., 2009, 8 (2), pp 974–981*

### Denator AB

Arvid Wallgrens Backe 20  
SE 413 46 Göteborg  
Sweden

Phone: +46 (0)31 412 890

Fax: +46 (0)31 412 840

inform@denator.com

**www.denator.com**

© 2010 Denator AB.

Stabilizor and Maintainor are trademarks of Denator AB. All third party trademarks are the property of their respective owners

User Manual, Denator 60007 REV 4

To be used together with DKT0001.